

1.10011.0001  
1.10011.0002MQuant™  
Peroxide Test

## 1. Method

Peroxidase transfers peroxide oxygen to an organic redox indicator. This produces a blue oxidation product. The peroxide concentration is measured **semi-quantitatively** by visual comparison of the reaction zone of the test strip with the fields of a color scale.

## 2. Measuring range and number of determinations

Measuring range / color-scale graduation	Number of determinations
0.5 - 2 - 5 - 10 - 25 mg/l H <sub>2</sub> O <sub>2</sub>	25 (Cat. No. 1.10011.0002) or 100 (Cat. No. 1.10011.0001)

## 3. Applications

This test measures inorganic peroxides in aqueous solutions and organic solvents. Polymeric peroxides are not at all or only incompletely measured.

**Sample material:**

Simple ethers  
UHT milk  
Pickling and copper-stripping baths  
Bleaching and oxidizing agents (paper and textile industries)  
Disinfectant and rinsing solutions (e. g. food technology, laundries)  
Swimming-pool water

## 4. Influence of foreign substances

This was checked in solutions with 12.5 and 0 mg/l H<sub>2</sub>O<sub>2</sub>. The determination is not yet interfered with up to the concentrations of foreign substances given in the table.

Concentrations of foreign substances in mg/l			
CrO <sub>4</sub> <sup>2-</sup>	10	IO <sub>3</sub> <sup>-</sup>	40
[Fe(CN) <sub>6</sub> ] <sup>4-</sup>	10	MnO <sub>4</sub> <sup>-</sup>	2
[Fe(CN) <sub>6</sub> ] <sup>3-</sup>	10	S <sub>2</sub> O <sub>8</sub> <sup>2-</sup>	20
Hg <sup>+</sup>	250	VO <sub>3</sub> <sup>-</sup>	5

## 5. Reagents and auxiliaries

The test strips are stable up to the date stated on the pack when stored closed at +2 to +8 °C.

**Failure to adhere to the storage temperature of +2 to +8°C will lower the shelf life of the test strips and the accuracy of the measuring values.**

**Package contents:**

Tube containing 25 test strips  
(Cat. No. 1.10011.0002)  
or  
containing 100 test strips  
(Cat. No. 1.10011.0001)

**Other reagents:**

MColorpHast™ Universal indicator strips pH 0 - 14, Cat. No. 109535  
Sodium acetate anhydrous GR for analysis, Cat. No. 106268  
Hydrochloric acid 1 mol/l Titripur®, Cat. No. 109057  
Diethyl ether for analysis EMSURE®, Cat. No. 100921  
Hydrogen peroxide 30 % H<sub>2</sub>O<sub>2</sub> (Perhydro®) GR for analysis, Cat. No. 107209

## 6. Preparation

- Samples containing more than 25 mg/l H<sub>2</sub>O<sub>2</sub> must be diluted with distilled water or peroxide-free ether.
- The pH of the aqueous sample must be within the range 2 - 12.  
If necessary, buffer the sample with sodium acetate or, respectively, adjust the pH with hydrochloric acid.

## 7. Procedure

**Protect the reaction zones from light (also during the reaction time)!**

**Determination in aqueous solutions:**

Immerse the reaction zone of the test strip in the pre-treated sample (15 - 30 °C) for 1 sec.

Allow excess liquid to run off via the long edge of the strip onto an absorbent paper towel and **after 15 sec** determine with which color field on the label the color of the reaction zone coincides most exactly.

Read off the corresponding result in mg/l H<sub>2</sub>O<sub>2</sub>.

**Notes on the measurement:**

- Every blue colouration **within 3 min** can be interpreted as a positive result.
- If the color of the reaction zone is equal to or more intense than the darkest color on the scale or if another color emerges, repeat the measurement using **fresh** samples diluted with distilled water or, respectively, peroxide-free ether until a value of less than 25 mg/l H<sub>2</sub>O<sub>2</sub> is obtained.

**Determination in organic solvents (readily volatile ethers):**

Immerse the reaction zone of the test strip in the pre-treated sample (15 - 30 °C) for 1 sec.

After the solvent has evaporated (gently fan the strip back and forth for 3 - 30 sec), humidify the reaction zone for 1 sec with 1 drop of distilled water and allow excess liquid to run off via the long edge of the strip onto an absorbent paper towel.

**After 15 sec** assess the color of the reaction zone.

**Note on the measurement:**

**It is recommended to treat the measurement results obtained in organic solvents only as guideline values, since the color in appearance and intensity may vary depending on the solvent medium. In this connection every blue coloration of the reaction zone indicates that peroxide is present.**

## 8. Method control

To check test strips and handling:

Make up 5.0 ml of Perhydro® (H<sub>2</sub>O<sub>2</sub> 30 % ± 333 000 mg/l H<sub>2</sub>O<sub>2</sub>) to 500 ml with distilled water and mix.

Take 1.5 ml of this solution, make up to 500 ml with distilled water, and mix. Subsequently analyze **immediately (solution is not stable)** as described in section 7. The content of H<sub>2</sub>O<sub>2</sub> determined should be 10 mg/l.

Additional notes see under [www.qa-test-kits.com](http://www.qa-test-kits.com).

## 9. Note

**Reclose the tube containing the test strips immediately after use.**

